

MS400C Hotplate Magnetic Stirrer




Features

- Ceramic coated top plate is easy to clean
- Large LCD display clearly indicates running time, temperature and stirring speed
- Push-button speed control prevents the stir bar from stalling accidentally
- Connected temperature probe automatically controls heating and maintains a constant temperature
- Single point temperature calibration



Specifications

Model	MS400C
Maximum Stirring Volume	2000 mL (H ₂ O)
Stirring Speed	0...2000 rpm
Timer Range	0...999 minutes
Hotplate Temperature	Max. 300°C
Top Plate Size	135(L)×135(W)mm
Top Plate Material	Ceramic coating
Display	Custom LCD (95 × 35 mm)
Power Supply	AC 220V/50Hz
Environmental Conditions	< 50°C, relative humidity < 80%
Dimensions and Weight	230(L)×180(W)×120(H)mm, 2.2 kg



Ordering Information

MS400C: Stirrer, stir bar, temperature probe, stand set